



Article available online at: <http://www.rt-image.com/productpreview2-7>

CR scanner is designed for full DICOM connectivity

## Product Preview

02.07.11



**The FireCR medical CR scanner is designed for full DICOM Connectivity. (Courtesy: 3DISC)**

CR manufacturer **3DISC** has introduced a line of upgradeable CR scanners. The **FireCR** medical CR scanner can be configured for most medical and chiropractic applications and is designed for full DICOM connectivity, which allows users to capture high-quality X-ray images.

The FireCR scanner can be upgraded to higher speeds with a smart card update to give users flexibility to start with an entry-level scanner and increase productivity. The FireCR medical product line includes a range of compact, portable CR readers for specific clinical applications.

**For more information, go to [www.3-disc.com](http://www.3-disc.com) (<http://www.3-disc.com>).**

### PORTABLE WORKSTATION

**Nextcomputing**, Nashua, N.H., a provider of portable and small footprint workstations and servers, has announced its newest portable workstation, the **Radius EX**. The Radius EX offers desktop-class computing performance in a transportable, all-in-one design. The workstation can run high-end software applications outside of typical office environments or in multiple locations.

Highlights of the Radius EX include: a compact, all-in-one chassis with a built-in 17-inch 1920x1200 display; durable aluminum alloy chassis with a padded leather handle; a single Intel Core i7 processor or dual Intel Xeon 5600 series processors; PCI Express or PCI expansions slots; 48GB of RAM; SATA or SAS 2.5-inch fixed internal hard drives, for

more than 7TB (terabytes) of high-speed storage, which is capable of more than 1GB/s (gigabyte-per-second) of data throughput.

**For more information, go to [www.nextcomputing.com](http://www.nextcomputing.com) (<http://www.nextcomputing.com>).**

## **EXECUTIVE DASHBOARD**

**Novarad**, American Fork, Utah, has introduced the executive dashboard **NovaDash**. An executive dashboard presents historical, current, and predictive views of operations to deliver insight into the health of a practice. NovaDash supports better decision-making by presenting data that tracks productivity, recognizes operational efficiencies, identifies and corrects trends before they become an issue, and monitors and verifies compliance.

The NovaRad dashboard refreshes with current data from diverse data sources and can interface with any HIS, and billing or compliant databases. The dashboard generates customized predefined data elements, such as stat reads, unread studies, financial data, or compliance reporting. This information enables a practice to take corrective measures immediately.

**For more information, go to [www.novarad.net](http://www.novarad.net) (<http://www.novarad.net>).**

---

Copyright © 2011, All Rights Reserved, Valley Forge Publishing Group