

Software/Hardware Requirements

Software Requirements

- Windows 10 (Excluding Windows 10 S, now defunct)
- > 2 GB of free disk space
- Open USB 3.0 port

Recommended Hardware

- 12 GB of RAM (DDR4 or better)
- ≥ 4 Core Intel i7 or AMD Ryzen Processor (at ≥ 2.8 GHz clock or ≥ 3.2 GHz turbo clock)
- NVIDIA Graphics Processor* - **5,000 GFLOPs** or greater processing power
 - 900 Series and below: **NOT** Compatible
 - 1000 Series:
 - 1070 or Greater
 - Mobile Series:
 - 1070 or Greater

* AMD GPU is *NOT* compatible with QuantorClinic and the Heron IOS.

Recommended High-Performance Hardware

- 16 GB of RAM (DDR4 or better)
- ≥ 6 Core Intel i7 or AMD Ryzen Processor (at ≥ 3.0 GHz clock or ≥ 3.6 GHz turbo clock)
- NVIDIA Graphics Processor* - **8,000 GFLOPs** or greater processing power
 - 900 Series and below: **NOT** Compatible
 - 1000 Series:
 - 1080 or Greater
 - Mobile Series:
 - 1080 or Greater

Examples of Compatible Laptop Computers

Brand:	Model #:
Alienware	AW17R5-7811
Alienware	AW17R5-7798
Alienware	AW17R5-9191
Alienware	AW17R5-9729