

Software/Hardware Requirements

Software Requirements

- Windows 10 (Excluding Windows 10 S, now defunct)
- ≥ 100 GB of free disk space
- 1x USB 3.0 port

Recommended Hardware

- 16 GB of RAM (DDR4 or better)
- ≥ 4 Core Intel i7 or greater (e.g. i7 7800) with ≥ 2.1 GHz base clock, ≥ 4.0 GHz Turbo/Boost Clock
- NVIDIA GeForce Graphics Processor
 - 10 Series (GTX): 1070 or greater - At least 6 GB video memory
 - 20 Series (RTX): 2060 or greater - At least 6 GB video memory

* AMD GPU is *NOT* compatible with QuantorClinic and the Heron IOS.

